Phyloseq tutorial

Daniel Vaulot 23 01 2018

Contents

| 1 | Aim | 1 | | | | | | | | | |
|---|-----|---|--|--|--|--|--|--|--|--|--|
| | 1.1 | Phyloseq R library | | | | | | | | | |
| | 1.2 | Data | | | | | | | | | |
| | 1.3 | References for data | | | | | | | | | |
| | 1.4 | To be added | | | | | | | | | |
| 2 | Pre | requisites to be installed | | | | | | | | | |
| 3 | | cript description | | | | | | | | | |
| | 3.1 | Load necessary libraries | | | | | | | | | |
| | 3.2 | Read the data and create phyloseq objects | | | | | | | | | |
| | 3.3 | Bar graphs | | | | | | | | | |
| | 3.4 | Heatmaps | | | | | | | | | |
| | 3.5 | Alpha diversity | | | | | | | | | |
| | 3.6 | Ordination | | | | | | | | | |
| | 3.7 | Network analysis | | | | | | | | | |


1 Aim

This document explains the use of the phyloseq R library to analyze metabarcoding data.

1.1 Phyloseq R library

- Phyloseq web site : $\label{eq:https:/joey711.github.io/phyloseq/index.html} \text{ }$
- $\bullet\,$ See in particular tutorials for
 - -importing data: https://joey711.github.io/phyloseq/import-data.html
 - -heat maps: https://joey711.github.io/phyloseq/plot_heatmap-examples.html

1.2 Data



This tutorial uses a reduced metabarcoding dataset obtained by C. Ribeiro and A. Lopes dos Santos. This dataset originates from the CARBOM cruise in 2013 off Brazil and corresponds to the 18S V4 region amplified on flow cytometry sorted samples (see pptx file for details) and sequenced on an Illumina run 2*250 bp analyzed with mothur.

1.3 References for data

- Gérikas Ribeiro, C., Lopes dos Santos, A., Marie, D., Helena Pellizari, V., Pereira Brandini, F., and Vaulot, D. (2016). Pico and nanoplankton abundance and carbon stocks along the Brazilian Bight. PeerJ 4, e2587. doi:10.7717/peerj.2587.
- Gérikas Ribeiro, C., Marie, D., Lopes dos Santos, A., Pereira Brandini, F., and Vaulot, D. (2016). Estimating microbial populations by flow cytometry: Comparison between instruments. Limnol. Oceanogr. Methods 14, 750–758. doi:10.1002/lom3.10135.
- Gérikas Ribeiro C, Lopes dos Santos A, Marie D, Brandini P, Vaulot D. (2018). Relationships between photosynthetic eukaryotes and nitrogen-fixing cyanobacteria off Brazil. ISME J in press.

1.4 To be added

• Exercices

2 Prerequisites to be installed

- R studio: https://www.rstudio.com/products/rstudio/download/#download
- Download and install the following libraries by running under R studio the following lines

```
install.packages("dplyr")  # To manipulate dataframes
install.packages("readxl")  # To read Excel files into R

install.packages("ggplot2")  # for high quality graphics
```

```
source("https://bioconductor.org/biocLite.R")
biocLite("phyloseq")
```

3 Script description

• Load the script file : $Phyloseq_tutorial.R$

3.1 Load necessary libraries

```
library("phyloseq")
library("ggplot2")
                        # graphics
library("readxl")
                        # necessary to import the data from Excel file
library("dplyr")
                        # filter and reformat data frames
## Attaching package: 'dplyr'
## The following objects are masked from 'package:stats':
##
##
       filter, lag
## The following objects are masked from 'package:base':
##
##
       intersect, setdiff, setequal, union
```

3.2 Read the data and create phyloseq objects

Three tables are needed

- OTU
- Taxonomy
- Samples

They are read from a single Excel file where each sheet contains one of the tables

OTU table

| 4 | A | В | С | D | E | F | G | H | 1 | J | K | L | M | N | 0 | P | Q | R | S | Т | U |
|----|--------|-------|------|-------|------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| | otu | | | | | | | | | | | | | | | | | | | | |
| 1 | | X10n | X10p | X11n | X11p | X120n | X120p | X121n | X121p | X122n | X122p | X125n | X125p | X126n | X126p | X127n | X13n | X13p | X140n | X140p | X141n |
| 2 | Otu001 | 13679 | 6292 | 42 | 2500 | 18850 | 5 | 43 | 7138 | 9432 | 10541 | 9 | 9772 | 1388 | 7 | 31538 | 38 | 2338 | 23 | 9 | 135 |
| 3 | Otu002 | 18 | 7134 | 38 | 9830 | 45 | 61420 | 182 | 23751 | 36 | 11 | 4535 | 3502 | 11018 | 5473 | 26 | 14411 | 38 | 19018 | 12 | 3080 |
| 4 | Otu003 | 9939 | 8983 | 31 | 13 | 24620 | 19 | 19 | 16 | 12502 | 3831 | 4621 | 2240 | 9924 | 4052 | 9292 | 18 | 0 | 37 | 7 | 3680 |
| 5 | Otu004 | 3675 | 4234 | 24 | 22 | 11 | 16 | 32967 | 35 | 6 | 18 | 6908 | 5 | 16 | 8702 | 24 | 11 | 37717 | 0 | 25 | 4196 |
| 6 | Otu005 | 0 | 5 | 0 | 7 | 0 | 8 | 0 | 16 | 20166 | 0 | 0 | 2 | 5 | 8 | 2 | 16 | 0 | 13 | 0 | (|
| 7 | Otu006 | 0 | 8 | 0 | 0 | 0 | 8 | 0 | 0 | 5 | 3 | 3 | 0 | 0 | 9 | 0 | 5 | 4 | 0 | 0 | |
| 8 | Otu007 | 4587 | 518 | 4 | 386 | 8775 | 5 | 6 | 1102 | 14336 | 0 | 0 | 3626 | 51 | 0 | 6 | 12 | 0 | 10 | 0 | 395 |
| 9 | Otu008 | 1 | 8 | 2 | 4408 | 3 | 29 | 6 | 12355 | 0 | 0 | 0 | 0 | 0 | 9 | 3 | 1588 | 0 | 6 | 3 | 3 |
| 10 | Otu009 | 115 | 914 | 3 | 325 | 0 | 629 | 1 | 834 | 5 | 0 | 1354 | 2108 | 1117 | 67 | 0 | 2010 | 1897 | 11227 | 1 | 3 |
| 11 | Otu010 | 780 | 8 | 23810 | 12 | 3279 | 0 | 12 | 7 | 3027 | 0 | 2 | 4156 | 0 | 0 | 18 | 0 | 0 | 0 | 0 | C |
| 12 | Otu011 | 0 | 3 | 2 | 2 | 0 | 13 | 5 | 5 | 4 | 7 | 3081 | 11 | 4 | 6804 | 0 | 3 | 11 | 0 | 5 | 0 |
| 13 | Otu012 | 0 | 0 | 0 | 6 | 0 | 0 | 0 | 16 | 3 | 0 | 0 | 0 | 0 | 0 | 0 | 17 | 0 | 6 | 0 | C |
| 14 | Otu013 | 6321 | 2471 | 2 | 0 | 12 | 3 | 0 | 0 | 4 | 20272 | 0 | 15 | 9 | 0 | 5 | 0 | 11 | 0 | 14 | (|
| 15 | Otu014 | 0 | 82 | 4 | 3304 | 1 | 1667 | 4 | 9233 | 13 | 3 | 0 | 2707 | 0 | 0 | 3 | 4806 | 9 | 3 | 5 | (|
| 16 | Otu015 | 0 | 12 | 0 | 3 | 7 | 25 | 1 | 6 | 10 | 0 | 4 | 2772 | 1 | 3 | 0 | 2 | 0 | 10 | 13 | 8052 |
| 17 | Otu016 | 1 | 0 | 0 | 9 | 5 | 0 | 0 | 14 | 0 | 0 | 0 | 0 | 2654 | 0 | 0 | 6 | 1 | 1 | 0 | C |
| 18 | Otu017 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 17 | 8 | 0 | 0 | 0 | 0 | 0 | 17 | 24 | 48 | 35210 | 4 |
| 19 | Otu018 | 1 | 0 | 9 | 911 | 0 | 0 | 15 | 2702 | 6 | 4 | 342 | 2217 | 606 | 0 | 13 | 3846 | 4 | 6 | 8513 | |
| 20 | Otu019 | 0 | 0 | 13 | 0 | 0 | 0 | 29 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 11 | 0 | 0 | 5 | 4 | C |
| 21 | Otu020 | 425 | 0 | 1 | 0 | 1706 | 0 | 8447 | 1 | 0 | 0 | 0 | 0 | 0 | 26 | 0 | 0 | 3490 | 0 | 2620 | (|
| 22 | Otu021 | 0 | 4 | 0 | 0 | 0 | 10 | 0 | 0 | 0 | 0 | 2 | 0 | 0 | 4 | 0 | 0 | 0 | 0 | 0 | - 4 |
| 23 | Otu022 | 0 | 0 | 0 | 4987 | 0 | 0 | 0 | 6 | 90 | 1 | 1 | 524 | 0 | 467 | 0 | 4 | 8 | 6198 | 0 | |
| 24 | Otu023 | 4 | 0 | 1 | 0 | 0 | 3 | 0 | 0 | 0 | 0 | 0 | 0 | 3351 | 3 | 0 | 3910 | 1 | 2 | 3 | (|
| 25 | Otu024 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 17 | 0 | 0 | 0 | 0 | 0 | 2 | 1 | (|
| 26 | Otu025 | 69 | 0 | 0 | 0 | 290 | 0 | 0 | 0 | 21 | 0 | 118 | 2 | 9 | 513 | 2 | 0 | 0 | 2 | 0 | (|
| 27 | Otu026 | 0 | 2 | 0 | 0 | 0 | 3 | 0 | 0 | 0 | 1 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | C |
| 28 | Otu027 | 6 | 2304 | 0 | 0 | 5 | 0 | 0 | 0 | 57 | 4 | 0 | 14529 | 9597 | 2 | 6 | 0 | 0 | 0 | 0 | C |

Taxon table

| | Α | В | С | D | E | F | G | н |
|---|--------|-----------|----------------|----------------|---------------------|-----------------------|-----------------------------|------------------------------|
| | otu | Domain | Supergroup | Division | Class | Order | Family | Genus |
| | Otu001 | Eukaryota | Archaeplastida | Chlorophyta | Mamiellophyceae | Mamiellales | Bathycoccaceae | Ostreococcus |
| | Otu002 | Eukaryota | Hacrobia | Haptophyta | Prymnesiophyceae | Prymnesiophyceae_X | Braarudosphaeraceae | UCYN_A1_host |
| | Otu003 | Eukaryota | Archaeplastida | Chlorophyta | Mamiellophyceae | Mamiellales | Bathycoccaceae | Bathycoccus |
| | Otu004 | Eukaryota | Alveolata | Dinophyta | Dinophyceae | Dinophyceae_X | Dinophyceae_X | Prorocentrum |
| | Otu005 | Eukaryota | Stramenopiles | Ochrophyta | Bacillariophyta | Mediophyceae | Mediophyceae_X | Thalassiosira |
| | Otu006 | Eukaryota | Stramenopiles | Ochrophyta | Bacillariophyta | Bacillariophyceae | Bacillariophyceae_X | Pseudo_nitzschia |
| | Otu007 | Eukaryota | Stramenopiles | Ochrophyta | Pelagophyceae | Pelagophyceae_X | Pelagophyceae_X | Pelagomonas |
| | Otu008 | Eukaryota | Alveolata | Dinophyta | Dinophyceae | Dinophyceae_X | Dinophyceae_X | Dinophyceae_X |
| | Otu009 | Eukaryota | Hacrobia | Haptophyta | Prymnesiophyceae | Prymnesiales | Chrysochromulinaceae | Chrysochromulina |
| | Otu010 | Eukaryota | Opisthokonta | Metazoa | Craniata | Craniata_X | Craniata_XX | Craniata_XX_unclassified |
| | Otu011 | Eukaryota | Stramenopiles | Ochrophyta | Chrysophyceae | Chrysophyceae_X | Chrysophyceae_Clade_C | Chrysophyceae_Clade_C_X |
| (| Otu012 | Eukaryota | Alveolata | Dinophyta | Dinophyceae | Dinophyceae_X | Dinophyceae_X | Gonyaulax |
| | Otu013 | Eukaryota | Alveolata | Dinophyta | Syndiniales | Syndiniales_Group_III | Syndiniales_Group_III_X | Syndiniales_Group_III_X |
| | Otu014 | Eukaryota | Stramenopiles | Ochrophyta | Chrysophyceae | Chrysophyceae_X | Chrysophyceae_Clade_G | Chrysophyceae_Clade_G_X |
| | Otu015 | Eukaryota | Alveolata | Dinophyta | Dinophyceae | Dinophyceae_X | Dinophyceae_X | Dinophyceae_X |
| 1 | Otu016 | Eukaryota | Hacrobia | Centroheliozoa | Centroheliozoa_X | Pterocystida | Pterocystida_X | Pterocystida_X |
| | Otu017 | Eukaryota | Opisthokonta | Fungi | Basidiomycota | Agaricomycotina | Agaricomycetes | Hyphodontia |
| | Otu018 | Eukaryota | Stramenopiles | Ochrophyta | Dictyochophyceae | Dictyochophyceae_X | Pedinellales | Pedinellales_X |
| | Otu019 | Eukaryota | Opisthokonta | Fungi | Basidiomycota | Agaricomycotina | Agaricomycetes | Itersonilia |
| 1 | Otu020 | Eukaryota | Hacrobia | Haptophyta | Prymnesiophyceae | Prymnesiophyceae_X | Braarudosphaeraceae | Braarudosphaera |
| | Otu021 | Eukaryota | Alveolata | Dinophyta | Dinophyceae | Dinophyceae_X | Dinophyceae_X | Dinophyceae_X |
| 1 | Otu022 | Eukaryota | Hacrobia | Haptophyta | Prymnesiophyceae | Prymnesiophyceae_X | Prymnesiophyceae_X | Syracosphaera |
| 1 | Otu023 | Eukaryota | Stramenopiles | Ochrophyta | Bacillariophyta | Bacillariophyceae | Bacillariophyceae_X | Bacillariophyceae_X |
| 1 | Otu024 | Eukaryota | Archaeplastida | Streptophyta | Klebsormidiophyceae | Klebsormidiophyceae_X | Klebsormidiophyceae_XX | Klebsormidium |
| 1 | Otu025 | Eukaryota | Archaeplastida | Chlorophyta | Mamiellophyceae | Mamiellales | Mamiellaceae | Micromonas |
| 1 | Otu026 | Eukaryota | Stramenopiles | Ochrophyta | Bacillariophyta | Bacillariophyceae | Bacillariophyceae_X | Cylindrotheca |
| | Otu027 | Eukaryota | Alveolata | Dinophyta | Dinophyceae | Suessiales | Suessiales_X | Karlodinium |
| | Otu028 | Eukaryota | Hacrobia | Haptophyta | Prymnesiophyceae | Isochrysidales | Noelaerhabdaceae | Emiliania |
| | Otu029 | Eukaryota | Opisthokonta | Fungi | Ascomycota | Saccharomycotina | Saccharomycetales | Debaryomyces |
| | Otu030 | Eukaryota | Hacrobia | Cryptophyta | Cryptophyceae | Cryptophyceae_X | Cryptomonadales | Teleaulax |
| | Otu031 | Eukaryota | Alveolata | Dinophyta | Syndiniales | Syndiniales_Group_I | Syndiniales_Group_I_Clade_1 | Syndiniales_Group_I_Clade_1_ |
| ı | Otu032 | Eukaryota | Archaeplastida | Chlorophyta | Prasino_Clade_VII | Prasino_Clade_VII_X | Prasino_Clade_VII_A | Prasino_Clade_VII_A_4_X |
| | | | | | | | | |

Sample table

| 1 | Α | В | С | D | E | F | G | н | 1 | J | K | L |
|----|--------|----------|-----------------|-----------|-----------|------------|---------------|----------|---------|-------|----------|-----------|
| 1 | sample | fraction | Select 18S nifH | total 18S | total 16S | total nifH | sample number | transect | station | depth | latitude | longitude |
| 2 | X10n | Nano | Yes | 53230 | 8772 | 36 | 10 | 1 | 81 | 140 | -27.42 | -44.72 |
| 3 | X10p | Pico | Yes | 47390 | 4448 | 6241 | 10 | 1 | 81 | 140 | -27.42 | -44.72 |
| 4 | X11n | Nano | No | 24007 | 6193 | 3772 | - 11 | 1 | 85 | 110 | -26.8 | -45.3 |
| 5 | X11p | Pico | Yes | 31899 | 14 | 10201 | 11 | 1 | 85 | 110 | -26.8 | -45.3 |
| 6 | X120n | Nano | Yes | 70455 | 5292 | 93 | 120 | 2 | 96 | 5 | -27.39 | -47.82 |
| 7 | X120p | Pico | Yes | 76182 | 53272 | 23147 | 120 | 2 | 96 | 5 | -27.39 | -47.82 |
| 8 | X121n | Nano | Yes | 52401 | 5958 | 26838 | 121 | 2 | 96 | 30 | -27.39 | -47.82 |
| 9 | X121p | Pico | Yes | 71785 | 10993 | 23706 | 121 | 2 | 96 | 30 | -27.39 | -47.82 |
| 10 | X122n | Nano | Yes | 78740 | 11730 | 15543 | 122 | 2 | 96 | 50 | -27.39 | -47.82 |
| 11 | X122p | Pico | Yes | 37364 | 11817 | 11045 | 122 | 2 | 96 | 50 | -27.39 | -47.82 |
| 12 | X125n | Nano | Yes | 27381 | 9 | 14331 | 125 | 2 | 98 | 5 | -27.59 | -47.39 |
| 13 | X125p | Pico | Yes | 55179 | 10419 | 21461 | 125 | 2 | 98 | 5 | -27.59 | -47.39 |
| 14 | X126n | Nano | Yes | 65714 | 15 | 16929 | 126 | 2 | 98 | 50 | -27.59 | -47.39 |
| 15 | X126p | Pico | Yes | 30406 | 3 | 10140 | 126 | 2 | 98 | 50 | -27.59 | -47.39 |
| 16 | X127n | Nano | Yes | 60610 | 9 | 11493 | 127 | 2 | 98 | 85 | -27.59 | -47.39 |
| 17 | X13n | Nano | Yes | 46001 | 33 | 21316 | 13 | 1 | 86 | 105 | -26.33 | -45.41 |
| 18 | X13p | Pico | Yes | 59626 | 7217 | 11954 | 13 | 1 | 86 | 105 | -26.33 | -45.41 |
| 19 | X140n | Nano | Yes | 48126 | 10428 | 25286 | 140 | 2 | 101 | 5 | -27.79 | -46.96 |
| 20 | X140p | Pico | Yes | 46569 | 10448 | 12301 | 140 | 2 | 101 | 5 | -27.79 | -46.96 |
| 21 | X141n | Nano | Yes | 30081 | 6394 | 21302 | 141 | 2 | 101 | 60 | -27.79 | -46.96 |
| 22 | X141p | Pico | Yes | 64221 | 11318 | 10428 | 141 | 2 | 101 | 60 | -27.79 | -46.96 |
| 23 | X142n | Nano | Yes | 85219 | 23243 | 11753 | 142 | 2 | 101 | 110 | -27.79 | -46.96 |
| 24 | X142p | Pico | Yes | 89797 | 9553 | 17156 | 142 | 2 | 101 | 110 | -27.79 | -46.96 |
| | X155n | Nano | Yes | 54162 | 8237 | 20674 | 155 | 2 | 106 | 5 | -28.12 | -46.17 |
| 26 | X155p | Pico | Yes | 50782 | 7384 | 66172 | 155 | 2 | 106 | 5 | -28.12 | -46.17 |
| 27 | X156n | Nano | Yes | 55065 | 11371 | 14447 | 156 | 2 | 106 | 60 | -28.12 | -46.17 |
| 28 | X156p | Pico | Yes | 43917 | 9665 | 16093 | 156 | 2 | 106 | 60 | -28.12 | -46.17 |
| 29 | X157n | Nano | Yes | 29078 | 4978 | 15532 | 157 | 2 | 106 | 100 | -28.12 | -46.17 |
| 30 | X157p | Pico | Yes | 51848 | 9139 | 15204 | 157 | 2 | 106 | 100 | -28.12 | -46.17 |
| 31 | X15n | Nano | Yes | 22468 | 2887 | 2678 | 15 | 1 | 87 | 105 | -26.22 | -45.48 |
| 32 | X15p | Pico | Yes | 78390 | 13813 | 1033 | 15 | 1 | 87 | 105 | -26.22 | -45.48 |
| 33 | X165n | Nano | Yes | 50732 | 15337 | 14706 | 165 | 2 | 114 | 5 | -28.65 | -44.99 |
| 34 | X165p | Pico | Yes | 48514 | 10902 | 39918 | 165 | 2 | 114 | 5 | -28.65 | -44.99 |
| 35 | X166n | Nano | Yes | 53412 | 3411 | 24442 | 166 | 2 | 114 | 60 | -28.65 | -44.99 |

```
otu_mat<- read_excel("CARBOM data.xlsx", sheet = "OTU matrix")
tax_mat<- read_excel("CARBOM data.xlsx", sheet = "Taxonomy table")
samples_df <- read_excel("CARBOM data.xlsx", sheet = "Samples")</pre>
```

```
Phyloseq objects need to have row.names
  • define the row names from the otu column
    row.names(otu_mat) <- otu_mat$otu</pre>
## Warning: Setting row names on a tibble is deprecated.
  • remove the column otu since it is now used as a row name
    otu_mat <- otu_mat %>% select (-otu)
  • Idem for the two other matrixes
 row.names(tax_mat) <- tax_mat$otu</pre>
## Warning: Setting row names on a tibble is deprecated.
  tax_mat <- tax_mat %>% select (-otu)
  row.names(samples_df) <- samples_df$sample</pre>
## Warning: Setting row names on a tibble is deprecated.
  samples_df <- samples_df %>% select (-sample)
Transform into matrixes otu and tax tables (sample table can be left as data frame)
  otu mat <- as.matrix(otu mat)</pre>
  tax_mat <- as.matrix(tax_mat)</pre>
Transform to phyloseg objects
  OTU = otu_table(otu_mat, taxa_are_rows = TRUE)
  TAX = tax_table(tax_mat)
  samples = sample_data(samples_df)
  carbom <- phyloseq(OTU, TAX, samples)</pre>
  carbom
## phyloseq-class experiment-level object
                                      [ 287 taxa and 55 samples ]
## otu table()
                  OTU Table:
## sample_data() Sample Data:
                                      [ 55 samples by 27 sample variables ]
                                      [ 287 taxa by 7 taxonomic ranks ]
## tax_table()
                  Taxonomy Table:
Visualize data
  sample_names(carbom)
  [1] "X10n"
                  "X10p"
                            "X11n"
                                     "X11p"
                                               "X120n"
                                                         "X120p"
                                                                  "X121n"
## [8] "X121p"
                  "X122n"
                            "X122p"
                                     "X125n"
                                               "X125p"
                                                         "X126n"
                                                                  "X126p"
## [15] "X127n"
                  "X13n"
                            "X13p"
                                     "X140n"
                                               "X140p"
                                                         "X141n"
                                                                  "X141p"
## [22] "X142n"
                  "X142p"
                            "X155n"
                                     "X155p"
                                               "X156n"
                                                         "X156p"
                                                                  "X157n"
## [29] "X157p"
                  "X15n"
                            "X15p"
                                     "X165n"
                                               "X165p"
                                                         "X166n"
                                                                  "X166p"
## [36] "X167n"
                  "X167p"
                            "X1n"
                                     "X1p"
                                               "X2n"
                                                         "X2p"
                                                                   "X3n"
                  "X5n"
                            "X5p"
                                                         "X9n"
                                                                   "X9p"
## [43] "X3p"
                                     "X7n"
                                               "X7p"
## [50] "tri01n" "tri01p" "tri02n" "tri02p" "tri03n" "tri03p"
```

"Class"

"Order"

"Supergroup" "Division"

"Genus"

rank_names(carbom)

[1] "Domain"

[6] "Family"


```
sample_variables(carbom)
  [1] "fraction"
                             "Select_18S_nifH"
                                                  "total 18S"
##
  [4] "total 16S"
                             "total nifH"
                                                  "sample_number"
##
## [7] "transect"
                             "station"
                                                  "depth"
## [10] "latitude"
                             "longitude"
                                                  "picoeuks"
## [13] "nanoeuks"
                             "bottom_depth"
                                                  "level"
                                                  "time"
## [16] "transect_distance" "date"
## [19] "phosphates"
                             "silicates"
                                                  "ammonia"
## [22] "nitrates"
                                                  "temperature"
                             "nitrites"
## [25] "fluorescence"
                             "salinity"
                                                  "sample_label"
Keep only samples to be analyzed
  carbom <- subset_samples(carbom, Select_18S_nifH =="Yes")</pre>
  carbom
## phyloseq-class experiment-level object
## otu_table()
                 OTU Table:
                                     [ 287 taxa and 54 samples ]
## sample_data() Sample Data:
                                     [ 54 samples by 27 sample variables ]
                                     [ 287 taxa by 7 taxonomic ranks ]
## tax_table()
                 Taxonomy Table:
Keep only photosynthetic taxa
  carbom <- subset_taxa(carbom, Division %in% c("Chlorophyta", "Dinophyta", "Cryptophyta",
                                                  "Haptophyta", "Ochrophyta", "Cercozoa"))
  carbom <- subset_taxa(carbom, !(Class %in% c("Syndiniales", "Sarcomonadea")))</pre>
  carbom
## phyloseq-class experiment-level object
## otu_table()
                 OTU Table:
                                     [ 205 taxa and 54 samples ]
                                     [ 54 samples by 27 sample variables ]
## sample_data() Sample Data:
## tax_table()
                 Taxonomy Table:
                                     [ 205 taxa by 7 taxonomic ranks ]
Normalize number of reads in each sample using median sequencing depth.
  total = median(sample_sums(carbom))
  standf = function(x, t=total) round(t * (x / sum(x)))
  carbom = transform_sample_counts(carbom, standf)
```

The number of reads used for normalization is 44903.

3.3 Bar graphs


Basic bar graph based on Division

plot_bar(carbom, fill = "Division")




Make the bargraph nicer by removing OTUs boundaries. This is done by adding ggplot2 modifier.

```
plot_bar(carbom, fill = "Division") +
geom_bar(aes(color=Division, fill=Division), stat="identity", position="stack")
```




Regroup together Pico vs Nano samples

```
carbom_fraction <- merge_samples(carbom, "fraction")
plot_bar(carbom_fraction, fill = "Division") +
geom_bar(aes(color=Division, fill=Division), stat="identity", position="stack")</pre>
```



Keep only Chlorophyta and use color according to genus. Do separate panels Pico vs Nano and Surface vs Deep samples.

```
carbom_chloro <- subset_taxa(carbom, Division %in% c("Chlorophyta"))
plot_bar(carbom_chloro, x="Genus", fill = "Genus", facet_grid = level~fraction) +
geom_bar(aes(color=Genus, fill=Genus), stat="identity", position="stack")</pre>
```




3.4 Heatmaps

A basic heatmap using the default parameters.

```
plot_heatmap(carbom, method = "NMDS", distance = "bray")
```


Warning: Transformation introduced infinite values in discrete y-axis




It is very very cluttered. It is better to only consider the most abundant OTUs for heatmaps. For example one can only take OTUs that represent at least 20% of reads in at least one sample. Remember we normalized all the sampples to median number of reads (total). We are left with only 33 OTUS which makes the reading much more easy.

```
carbom_abund <- filter_taxa(carbom, function(x) sum(x > total*0.20) > 0, TRUE)
  carbom abund
## phyloseq-class experiment-level object
## otu table()
                  OTU Table:
                                      [ 33 taxa and 54 samples ]
## sample_data() Sample Data:
                                      [ 54 samples by 27 sample variables ]
## tax_table()
                  Taxonomy Table:
                                      [ 33 taxa by 7 taxonomic ranks ]
  otu_table(carbom_abund)[1:8, 1:5]
## OTU Table:
                        [8 taxa and 5 samples]
##
                         taxa are rows
##
           X10n
                 X10p
                        X11p X120n X120p
## Otu001 13339
                 7346
                        3804 12662
                 8329 14958
## Otu002
             18
                                30 36206
## Otu003
           9692 10488
                          20 16537
                                       11
           3584
## Otu004
                  4943
                          33
                                 7
                                        9
## Otu005
              0
                     6
                          11
                                 0
                                        5
## Otu006
              0
                     9
                           0
                                 0
                                        5
## Otu007
                         587
                              5894
                                        3
           4473
                  605
## Otu008
                     9
                        6707
                                       17
  plot_heatmap(carbom_abund, method = "NMDS", distance = "bray")
```

Warning: Transformation introduced infinite values in discrete y-axis




It is possible to use different distances and different multivaraite methods. For example Jaccard distance and MDS and label OTUs with Class, order by Class. We can also change the Palette (the default palette is a bit ugly...).



Another strategy is to do a heatmap for a specific taxonomy group. For example we can taget the Chlorophyta and then label the OTUs using the Genus.

Warning: Transformation introduced infinite values in discrete y-axis




3.5 Alpha diversity

Plot Chao1 richness estimator and Shannon diversity estimator.

```
plot_richness(carbom, measures=c("Chao1", "Shannon"))
```


 $\mbox{\tt \#\#}$ Warning: Removed 54 rows containing missing values (geom_errorbar).



Regroup together samples from the same fraction.

```
plot_richness(carbom, measures=c("Chao1", "Shannon"), x="level", color="fraction")
```

Warning: Removed 54 rows containing missing values (geom_errorbar).



3.6 Ordination

Do multivariate analysis based on Bray-Curtis distance and NMDS ordination.


```
carbom.ord <- ordinate(carbom, "NMDS", "bray")</pre>
## Square root transformation
## Wisconsin double standardization
## Run 0 stress 0.2317058
## Run 1 stress 0.2488414
## Run 2 stress 0.2534458
## Run 3 stress 0.2359226
## Run 4 stress 0.2576921
## Run 5 stress 0.2576468
## Run 6 stress 0.2472519
## Run 7 stress 0.2510679
## Run 8 stress 0.2519793
## Run 9 stress 0.2302728
## ... New best solution
## ... Procrustes: rmse 0.1066564 max resid 0.3796017
## Run 10 stress 0.2450473
## Run 11 stress 0.2444403
## Run 12 stress 0.2348005
## Run 13 stress 0.2422467
## Run 14 stress 0.2566082
## Run 15 stress 0.2407624
## Run 16 stress 0.2542353
## Run 17 stress 0.2632189
## Run 18 stress 0.2585642
## Run 19 stress 0.2532353
## Run 20 stress 0.2449511
```

*** No convergence -- monoMDS stopping criteria:

20: stress ratio > sratmax


##

Plot \mathbf{OTUs}




A bit confusing, so make it more easy to visualize by breaking according to taxonomic division.

OTUs




Now display **samples** and enlarge the points to make it more easy to read.




Display both samples and OTUs but in 2 different panels.

biplot



3.7 Network analysis

Simple network analysis



This is quite confusing. Let us make it more simple by using only major OTUs

